

# HEINE BETA X Ophthalmoscope

Our most advanced direct ophthalmoscope in our 80-year history with AOS+ and the new focusCONTROL.

## Optimised for hygienic reprocessing

- Scratch-resistant window made of tempered glass
- smoothSURFACE and recessed USB-C connection

## Durable, lightweight and ergonomic design

- Solid, integrated aluminium frame
- Housing is thermo-moulded out of a special high-impact material blend optimised for hygienic reprocessing
- Precisely balanced for excellent ergonomics and secure handling

## Smart power source options

- Charge status indicator
- Automatic switch-off function
- Rechargeable with USB-C or charging station

## Superior refractive error compensation

Select the precise focus adjustment in 74 single dioptre steps from -36D to +38D

## Microprocessor controlled smart illumination control wheel

- Simple, one-thumb operation
- Turn wheel quickly for large changes in illumination intensity
- Turn wheel slowly for precise fine tuning of illumination intensity
- Continuous dimmability between 1 and 100%
- Push wheel for on|off
- Always turns on to the last used illumination intensity setting



# HEINE BETA X Ophthalmoscope

## Crisp, clear views of the fundus with the new AOS<sup>+</sup> optical system and LED<sup>HQ</sup>

- The AOS<sup>+</sup> Aspherical Optical System is even more effective in reducing corneal reflections through the optimised separation of illumination and observation paths
- LED<sup>HQ</sup> illumination for homogeneous, bright illumination with true colour rendering for high contrast and extremely detailed views
- Optical system is sealed to guard against dust

## The new swivel-in polarisation filter

Depending on requirements and the purpose of the examination:

- **With polarisation:** Reduces reflections in the cornea and stray light
- **Without polarisation:** Wall reflex (macular reflex) and retinal arteriolar light reflex (copper & silver wiring) clearly visible




Product data



## 7 different apertures

All apertures can also be used with red-free and polarisation filters



## The new HEINE focusCONTROL

Easily switch focus between the anterior segment of the eye and the fundus



HEINE Optotechnik GmbH & Co. KG

Dornierstr. 6, 82205 Gilching, Germany  
www.heine.com



We reserve the right to change specification without notice.

03|25. A-000.00.274 EN

QUALITY  
MADE IN GERMANY

