


Our best imaging, reimagined.



Enhance your HEINE dermatoscope with the SLR | SLM Camera Adapter. The combination of the HEINE DELTA 30 or HEINE DELTA 30 PRO, the HEINE SLR | SLM Camera Adapter and a DSLR or DSLM camera creates a comprehensive system for documenting consistently high-quality images.


Upgrade your high-end dermatoscope for even better imaging. The perfect solution, not only for scientific applications and the daily work of clinical photographers.

- ▶ Universally applicable for DELTA 30 PRO and DELTA 30
- ▶ Improved visibility of the smallest structures through digital post-magnification
- ▶ Highest optical resolution
- ▶ Best visibility of details
- ▶ Easy refocusing

Finding the right adapter:

The HEINE SLR | SLM camera adapter is available either with a bayonet connection for digital single-lens reflex (DSLR) cameras and digital single-lens mirrorless (DSLM) cameras or with a T2 universal connection.

When purchasing the SLR | SLM Camera Adapter with a T2 universal connection, a T2 adapter with bayonet camera connection is required to connect the digital camera to the SLR | SLM Camera Adapter.



These HEINE camera adapters including the bayonet adapters are available for the following digital cameras:

Brand	Bayonet adapter	HEINE SLR SLM Camera Adapter
Canon	DSLR EOS EF	K-000.34.141
	DSLM EOS R	K-000.34.142
	DSLM EOS M	K-000.34.143
NIKON	DSLR NIKON F	K-000.34.144
	DSLM NIKON Z	K-000.34.145
SONY	DSLM SONY E	K-000.34.146
Universal	Already available	K-000.34.140